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HDE *EVOLUTION*

Horse Dental Equipment

www.horse-dental-equipment.com

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Rechargeable lithium-ion battery

MP 144350

High performance
Medium Prismatic cell

Saft always supplies MP cells in assemblies or as customized battery system constructions

Benefits

- Extended autonomy and life for mobile systems
- A broad operating temperature range
- Recommended for ruggedized designs
- Easy integration into compact and light systems
- Very high energy density (344 Wh/l and 143 Wh/kg)
- Unrivalled low temperature performance

Key features

- Excellent charge recovery after long storage, even at high temperature
- Maintenance-free
- Non-restricted for transport/ Non-assigned to Class 9 according to the UN Recommendations on the transport of dangerous goods - Model Regulations
- Underwriters Laboratories Component Recognition (File Number MH 12609) (UL)

Main applications

- Mobile asset tracking
- Small UPS
- Soldier of the future equipment
- Portable radios
- Professional portable lighting
- Bar code readers
- Portable payment terminal

Electrical characteristics

Nominal voltage (0.5 A rate at 20°C)	3.75 V
Typical capacity (at 0.5 A 20°C 2.5 V cut-off)	2.6 Ah
Nominal energy	9.75 Wh

Mechanical characteristics (Sleeved 100 % charged cell)

Thickness <i>(Thickness tends to increase with cycling, typically obtained after 500 cycles. Consult Saft) (At beginning of life 14.6 mm)</i>	14.9 mm
Width (max)	43.9 mm
Height (max)	54.5 mm
Typical weight	68 g
Lithium equivalent content	0.78 g
Volume	28 cm ³

Operating conditions

Charge method	Constant Current/Constant Voltage
Charge voltage	4.20 +/- 0.05 V
Maximum recommended charge current* *	2.6 A (C rate)
Charge temperature range*	-20°C to +60°C
Charge time at 20°C	To be set as a function of the charge current: C rate ➡ 2 to 3 h C/2 rate ➡ 3 to 4 h C/5 rate ➡ 6 to 7 h
Maximum continuous discharge current**	5.0 A (~2C rate)
Pulse discharge current	up to 10 A (~4C rate)
Discharge cut-off voltage	2.5 V
Discharge temperature range*	-50°C to +60°C

* For optimized charging below 0°C, 60°C and discharging at -50°C, consult Saft.

** Electronic protection circuits within battery packs may limit the maximum charge/discharge current allowable. Consult Saft.


MP 144350

Battery assembly

In order to operate properly, individual Li-Ion cells are mechanically and electrically integrated in battery assemblies specific to each application. The battery assembly incorporates electronics for performance, thermal and safety management.

Technology

- Graphite-based anode
- Lithium Cobalt oxide-based cathode
- Electrolyte: organic solvents
- Built-in redundant safety protections (*shutdown separator, circuit breaker, safety vent*)
- Batteries assembled from MP cells feature an electronic protection circuit

Built-in protection devices ensure safety in case of:

- Exposure to heat
- Exposure to direct sunlight for extended periods of time
- Short circuit
- Overcharge
- Overdischarge

When handling Saft MP batteries:

- Do not solder directly to cell terminal
- Do not disassemble
- Do not remove the protection circuit
- Do not incinerate

Transportation and storage:

- Store in a dry place at a temperature preferably not exceeding 30°C
- For long-term storage, keep the battery within a (30 ± 15) % state of charge

**Saft
Specialty Battery Group**

www.saftbatteries.com

Doc. N° 54044-2-1009

Information in this document is subject to change without notice and becomes contractual only after written confirmation by Saft.


Published by the Communications Department.

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
Produced by Arthur Associates Limited.


Capacity versus current at + 20°C

Typical discharge profiles (0.5 A - C/5 rate)

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